


UTERINE ELEVAT **O.R.** PRO™




The smart way to elevate the uterus


The **Uterine ElevatOR PRO™** is an easy-to-use uterine manipulator that exposes the patient's pelvic anatomy for better definition and optimal access to surgical landmarks during laparoscopic and robotic hysterectomies.

Its curved design aligns to the angle of the sacral curve allowing for gentle insertion and effective manipulation of the uterus. When applying cephalad pressure, the Uterine ElevatOR PRO™ delineates the vaginal fornices and clearly displaces the ureters, creating a safe operative margin for the colpotomy incision. Pneumoperitoneum is maintained throughout the procedure.

The Uterine ElevatOR PRO™ is available in four cervical cups sizes to accommodate varied anatomy.



Intrauterine Balloon



Ergonomic Handle

CONTACT US

Asia Pacific
1/32 Silkwood Rise
Carrum Downs VIC
3201 Australia

T +61 3 9413 5555
F +61 3 9413 5556

North America
1625 South Tacoma Way
Tacoma WA 98409
USA

T +1 253 441 6534
F +1 866 333 9524

info@theorcompany.com
www.theorcompany.com

Uterine Manipulation

Simple, safe, effective

- Four cervical cup sizes (32mm, 35mm, 37mm, 40mm) allow for optimal setting of the cervix and identification of the vaginal fornices, thus accommodating different gynecological procedures and patient characteristics.
- Cervical cups, made out of Ethylene Tetrafluoroethylene (EFTE), offer high temperature resistance and strength. They can be used with any electrosurgical and harmonic energy sources.
- Intrauterine balloon (10cc) stabilizes the Uterine Elevator PRO™ within the uterine cavity, providing a firm hold and control of the uterus
- Anatomical (s-curved) design contours to the patient's anatomy and offers optimal range of motion of the uterus (anteversion, retroversion and lateral movement)
- Double-cup design of the Uterine Elevator PRO™ secures the cervical cup and occludes the vagina, preventing the loss of pneumoperitoneum throughout the procedure
- Ergonomic handle offers full control of the position of the uterus thus giving surgeons better access to key surgical landmarks during laparoscopic and robotic hysterectomies
- Dye injection port supports diagnostic laparoscopic, minilaparotomy and fertility exams
- Single-piece device, fully assembled, reduces procedural setup time and eliminates additional costs (extra components, repair sterilization)
- Not made with natural rubber latex



Order Information

REF/Order #	Description	UOM
UE-TVPRO-32	Uterine Elevator PRO™ – Small, 32mm/d cervical cup	Box of 8, sterile
UE-TVPRO-35	Uterine Elevator PRO™ – Medium, 35mm/d cervical cup	Box of 8, sterile
UE-TVPRO-37	Uterine Elevator PRO™ – Large, 37mm/d cervical cup	Box of 8, sterile
UE-TVPRO-40	Uterine Elevator PRO™ – XLarge, 40mm/d cervical cup	Box of 8, sterile

ABOUT US

The O.R. Company develops, supplies and markets high quality, innovative surgical devices - from niche products to proprietary surgical instruments - for minimally invasive and open surgery.

Through collaboration with surgeons, clinical technicians and nurses around the world, we identify surgical needs and advance them into effective and safe products such as the Uterine Elevator PRO™.

CONTACT US

Asia Pacific
1/32 Silkwood Rise
Carrum Downs VIC
3201 Australia
T +61 3 9413 5555
F +61 3 9413 5556

North America
1625 South Tacoma Way
Tacoma WA 98409
USA

T +1 253 441 6534
F +1 866 333 9524
info@theorcompany.com
www.theorcompany.com