

NEW - AlphaDigiCam Complete

- Optional HDMI output enables a digital output resolution up to 1080i for external monitors.
- HDMI enables the connection to a variety of monitors.
- USB archiving capabilities.
- LED light technology – long lifetime and reduced service cost.
- PS2 Keyboard connection.
- Intuitive operation.
- AlphaNovoCam II technology. Remote function via the included camera head.
- Use of standard scopes and light cables.



Camera head:

Aspect ratio: 4:3

- Image sensor: 1/3" CCD
- Immersible (using the attached cap)

System:

- Power supply: 100/240 VAC, 50-60 Hz
- Fuses: 2 Fine fuses, 5x20 mm, delay acting 3,15 AH (I=1500A)
- Power consumption: 90 VA
- Video resolution: 470 lines
- Sensitivity: < 1 lux (1/1.4)

Video system · PAL/NTSC

Monitor · 10,4" TFT, 800 x 600 dots resolution

Shutter range · PAL 1/50 - 1/100.000 sec; NTSC 1/60 - 1/100.000 sec

Lamp · LED lamp, comparable 50 W, 5500 K

Light guide connection · Gimmi, Storz

Video out terminal

· 2x Video; 2x S-Video; 1x Y-Pb-Pr (Component); 1x SDI; 1x USB 2.0; HDMI

Standard

Protection class: I

Application part: Type BF

IEC 601-1, EN 60601-1; classification acc. 93/42/EEC: I

Features

White balance, window functions, boost function, anti-moiree-filter

Freeze function

Documentation system

Dimensions · Control unit: 295 x 128 x 278 mm (W x H x D); 5 kg

Delivery contains · CCU, mains cable, user manual, , Y/C cable, BNC cable, USB Stick

NORMEDI NORDIC

Normedi Norge AS | Normedi AB | Normedi ApS | Oy Benite Ab

Ørsnesalléen 31

3120 Nøtterøy

Norway