


LIQUID SCOPE WARMOR™




The smart way to minimize scope fogging

The **Liquid Scope WarmOR™** is a single-use device that warms laparoscopes and endoscopes to body temperature. Activated in seconds it prevents condensation that causes the scope to fog during procedures.

Without the need for capital equipment or electrical connections, the **Liquid Scope WarmOR™** represents a cost effective and environment friendly scope warming and defogging solution. Available in two sizes, it accommodates scopes of all sizes.



Not activated



Click to activate

CONTACT US

Asia Pacific
122 Balmain Street
Richmond, VIC 3121
Australia

T +61 3 9413 5555
F +61 3 9413 5556


North America
1625 South Tacoma Way
Tacoma WA 98409
USA

T +1 253 441 6534
F +1 866 333 9524


info@theorcompany.com
www.theorcompany.com

Simple, instant heat source

- Affordable, single-use device
- Intuitive, click activation
- Quick to warm, reaches peak temperature (50°C +/- 5°C) in 5 minutes
- Keeps the heat for 40 minutes up to one hour
- Accommodates all sizes and lengths of scopes, used in different procedures
- Self-contained device that requires no batteries or electrical connections
- Eliminates fire, burns and hot water hazards and thus improves OR safety
- Sterile, natural rubber latex, DEHP and alcohol free
- Available in two sizes



Without - Fogged view



With - Clear view through a pre-warmed scope

Ordering Information

REF/Order #	Description	UOM
PS-500-M	Liquid Scope WarmOR™ – Short version (two buttons)	Box of 10, sterile
PS-500	Liquid Scope WarmOR™ – Long version (three buttons)	Box of 10, sterile

CONTACT US

Asia Pacific
122 Balmain Street
Richmond, VIC 3121
Australia

T +61 3 9413 5555
F +61 3 9413 5556

North America
1625 South Tacoma Way
Tacoma WA 98409
USA

T +1 253 441 6534
F +1 866 333 9524

info@theorcompany.com
www.theorcompany.com

ABOUT US

The O.R. Company develops, supplies and markets high quality, innovative surgical devices - from niche products to proprietary surgical instruments - for minimally invasive and open surgery.

Through collaboration with surgeons, clinical technicians and nurses around the world, we identify surgical needs and advance them into effective and safe products such as the Liquid Scope WarmOR™.