

The Smart Way to Prevent Leakage

ONE SIZE FITS ALL

SINGLE USE SEAL FOR ENDOSCOPY



The Endoscopic Seal is used to prevent fluid leakage and/or air aspiration through the working channel of a hysteroscope, cystoscope or ureteroscope during introduction and use of an instrument. The Endoscopic Seal features a funnel guided top for easy insertion of instruments, self-sealing technology for leak-free use with multiple instruments, and a luer-lock adapter compatible with any female luer-lock fitting.

ART.NR	DESCRIPTION	UNITS/BOX
351-905-105	Backflow Seal	10

Send your order to support@normedi.com