

The Smart Way to Minimize Fogging

ONE SIZE FITS ALL

LIQUID SCOPE WARMOR™



The **Liquid Scope WarmOR™** is a single-use device that warms laparoscopes and endoscopes to body temperature. Activated in seconds it prevents condensation that causes the scope to fog during procedures.

The Liquid Scope WarmOR eliminates the need for capital equipment or electrical connections resulting in a simple, cost-effective scope warming and fog prevention solution.

ART.NR	DESCRIPTION	UNITS/BOX
PS-500-M	Liquid Scope warmer	10

Send your order to support@normedi.com